


### FINOFIBRATE IMPURITY A

<b>Chemical name:</b> 4-Chloro-4'-hydroxybenzophenone	<b>CAS No. :</b> [42019-78-3]  <b>Structure:</b>  <p>The chemical structure shows two benzene rings connected by a carbonyl group (C=O). The left benzene ring has a hydroxyl group (-OH) at the para position. The right benzene ring has a chlorine atom (-Cl) at the para position.</p>
<b>Category :</b> Impurities / Metabolite	
<b>Molecular Formula :</b> C <sub>13</sub> H <sub>9</sub> ClO <sub>2</sub>	
<b>Molecular Weight :</b> 232.67	
<b>Product Code :</b> RBPL-F0401	