


## LANSOPRAZOLE IMPURITY

<b>Chemical name:</b> Lansoprazole sulphone N-oxide	
<b>Category :</b> Impurities / Metabolite	<b>CAS No. :</b> [953787-54-7]
<b>Molecular Formula :</b> C <sub>16</sub> H <sub>14</sub> F <sub>3</sub> N <sub>3</sub> O <sub>4</sub> S	<b>Structure:</b> 
<b>Molecular Weight :</b> 401.36	
<b>Product Code :</b> RBPL-L0402	